

Zeiss Elyra 7 with Lattice SIM

Location:

CELLIM, building C2, room 1.13a

Booking alias:

Elyra7-A2


Contact:

cellim@ceitec.muni.cz

jakub.pospisil@ceitec.muni.cz

petra.machackova@ceitec.muni.cz

T: 54949 7085/6133



Manual:



Reservations:

<https://booking.ceitec.cz/PlanningBoard.html>

Overview:

The Elyra 7 is a motorized inverted microscope with resolution beyond the Abbe diffraction limit of the conventional light microscope. The microscope is equipped with a super-resolution system that uses a specialized lattice illumination pattern for 3D structured illumination (SIM) images and single-molecule localization methods (SMLM) like STORM and PALM. SIM is a technology that can be used with standard fluorophores, making it an easily accessible choice for super-resolution results. The lattice illumination allows for faster, gentler, and deeper imaging of live samples, and can reach up to 255 FPS during time-lapse acquisition. The unit has an incubation chamber for live imaging with controlled CO₂ and temperature. For fast and precise Z-movement there is a piezo Z stage insert. The system yields up to 8-fold improvement in voxel resolution in the X, Y, and Z dimensions. The microscope is equipped with four solid-state lasers (405, 488, 561 and 642 nm). The signal detection is assured by two sCMOS cameras which enable simultaneous detection of two different fluorophores.

Specifications:

Objectives

Plan-Neofluoar 10x / 0.30 AIR,

C-Apochromat 40x / 1.20 W,

Plan-Apochromat 40x / 1.40 OIL,

Plan-Apochromat 63x / 1.46 OIL,

Plan-Apochromat 100x / 1.46 OIL,

Plan-Apochromat 100x / 1.57 OIL-HI

Lasers

405 nm, 488 nm, 561 nm, 642 nm

Lasers shared between Lattice SIM and SMLM

Camera

PCO edge sCMOS camera, 1280 × 1280;
 pixel size 6.5 μm × 6.5 μm

Software

Zen Black

Filters for fluorescence

F-set 77 HE – DAPI, GFP, mRFP, Alexa 633

Duolink filter sets

SR SOLO, SR DUO, SR Quad

Filters for SIM and SMLM

FLEX, FLEX min, PURE, QUAD only

Detailed description:

Duolink filter sets


Filter sets for SIM and SMLM



Examples:

Apotome visualization of Actin filaments and focal adhesion in Mouse Embryonic Fibroblast.

Cyan DAPI: DNA, Red Alexa Fluor 568, White Phalloidin
 Alexa 488: Actin filaments.
 Objective: Plan-Apochromat 40x / 1.4 Oil. Final image was obtained using ZEN SIM processing.



Mitochondria of endothelial cells stained with Mitotracker-RED.

Picture *A* was acquired using standard widefield microscopy. The same position was acquired using lattice SIM modality of Elyra 7 (picture *B*). There are clear resolution improvements and inner structures of mitochondria are revealed, otherwise invisible under a standard widefield microscope.

